

FIRM LOCK PLATES

Firm-lock metaphyseal plate 3.5 for distal, medial humerus

Firm-lock placas 3.5 para metáfisis, para húmero distal-medial

Firm-lock plaques 3.5 pour métaphyse, pour humérus distal médial



Stainless Steel	Size	Length mm	Titanium
Item Code			Item Code
651.006	6 Holes	67	651T.006
651.007	7 Holes	80	651T.007
651.008	8 Holes	93	651T.008
651.009	9 Holes	106	651T.009
651.010	10 Holes	119	651T.010
651.011	11 Holes	132	651T.011
651.012	12 Holes	145	651T.012
651.013	13 Holes	158	651T.013
651.014	14 Holes	171	651T.014
651.015	15 Holes	184	651T.015
651.016	16 Holes	197	651T.016

With LC undercuts. All holes counted.
 Sizes not listed above can also be manufactured on request.