

FIRM LOCK PLATES

Firm-lock reconstruction plate 4.5/5.0
 Firm-lock placas de reconstrucción 4.5/5.0
 Firm-lock plaques de reconstruction 4.5/5.0



Stainless Steel	Size	Length mm	Titanium
Item Code			Item Code
662.003	3 Holes	56	662T.003
662.004	4 Holes	75	662T.004
662.005	5 Holes	94	662T.005
662.006	6 Holes	113	662T.006
662.007	7 Holes	132	662T.007
662.008	8 Holes	151	662T.008
662.009	9 Holes	170	662T.009
662.010	10 Holes	189	662T.010
662.011	11 Holes	208	662T.011
662.012	12 Holes	227	662T.012
662.013	13 Holes	246	662T.013
662.014	14 Holes	265	662T.014
662.015	15 Holes	284	662T.015
662.016	16 Holes	303	662T.016
662.017	17 Holes	322	662T.017
662.018	18 Holes	341	662T.018
662.019	19 Holes	360	662T.019
662.020	20 Holes	379	662T.020

Sizes not listed above can also be manufactured on request.