

FIRM LOCK PLATES

Firm-lock proximal humerus plate 3.5
 Firm-lock placa para húmero proximal 3.5
 Firm-lock plaques 3.5 pour humérus proximal



Stainless Steel	Size	Length mm	Titanium
Item Code			Item Code
673.003	3 Holes	90	673T.003
673.004	4 Holes	102	673T.004
673.005	5 Holes	114	673T.005
673.006	6 Holes	126	673T.006
673.007	7 Holes	138	673T.007
673.008	8 Holes	150	673T.008
673.009	9 Holes	162	673T.009
673.010	10 Holes	174	673T.010
673.011	11 Holes	186	673T.011
673.012	12 Holes	198	673T.012

Proximal locking holes accept 3.5 mm Locking Screws.

Combi holes in shaft accept 3.5 mm Locking, 3.5 mm and 4.0 mm Cortex, and 4.0 mm Cancellous Bone Screws

Sizes not listed above can also be manufactured on request.