

FIRM LOCK PLATES

Firm-lock proximal diaphyseal humeral plate 3.5, long
 Firm-lock larga, placa diafisaria de húmero proximal 3.5
 Firm-lock plaques 3.5 longues pour humérus proximal diaphysaire



Stainless Steel	Size	Length mm	Titanium
Item Code			Item Code
674.003	3 Holes	106	674T.003
674.004	4 Holes	124	674T.004
674.005	5 Holes	142	674T.005
674.006	6 Holes	160	674T.006
674.007	7 Holes	178	674T.007
674.008	8 Holes	196	674T.008
674.009	9 Holes	214	674T.009
674.010	10 Holes	232	674T.010
674.011	11 Holes	250	674T.011
674.012	12 Holes	268	674T.012
674.013	13 Holes	286	674T.013
674.014	14 Holes	304	674T.014

Proximal locking holes accept 3.5 mm Locking Screws.

Combi holes in shaft accept 3.5 mm Locking, 3.5 mm and 4.0 mm Cortex, and 4.0 mm Cancellous Bone Screws

Sizes not listed above can also be manufactured on request.