


FIRM LOCK PLATES

Firm-lock dorsal distal radius T-plate 2.4, straight

Firm-lock placas 2.4 en T para radio distal, recta

Firm-lock plaques 2.4 en T pour radius distal, droite



Stainless Steel	Head	Shaft	Length	Titanium
Item Code				Item Code
90313.477	3 Holes	3 Holes	40	90313T.477
90313.478	3 Holes	4 Holes	49	90313T.478

*With LC undercuts. All holes in shaft counted.
 Sizes not listed above can also be manufactured on request.*